


# BLutein™ Dyes

A complete ophthalmic dye range designed to help free you from safety profile concerns



# Safer dyes are needed for **ophthalmic surgery**

- Vital dyes are commonly used to **improve visualisation** of delicate intraocular tissues during ophthalmic surgery<sup>1,2</sup>
- Although many dyes are currently available, they are **synthetic** and associated with **concerns about retinal toxicity**<sup>1,3,4</sup>
- To allow you to focus on offering the best possible patient outcomes, the **ideal dye** for eye surgery should offer:

1

**High safety profile for intraocular use**

2


**Reliable staining**

3

**Selective staining**

4

**Rapid elimination during removal**



# Introducing BLutein™ Dyes: The complete ophthalmic dye range designed to help free you from safety profile concerns

The BLutein™ range consists of dyes for both posterior segment staining and anterior capsule staining, so you always have a dye that meets your surgical needs



for staining  
the vitreous



for staining  
the ILM




for staining  
the ILM and ERM



for staining  
the anterior capsule

ERM: Epiretinal membrane; ILM: Internal limiting membrane.




# BLutein™ Dyes: The first and only dyes for eye surgery to contain lutein

The BLutein™ Dyes all contain **naturally sourced lutein**

- Lutein is a pigment **physiologically present in the human retina**, known to be important for overall eye health<sup>5-8</sup>

Lutein may also help **protect the retina** during surgery by:

- Absorbing damaging blue light before it reaches retinal cells through its blue light filter action<sup>9,10</sup>
- Acting as an antioxidant to neutralise free radicals that may be harmful to retinal cells<sup>11-13</sup>
- Reducing iatrogenic stress to the retinal tissue allowing faster functional recovery after surgery<sup>14</sup>




# The BLutein™ Dyes for ILM and ERM staining also contain highly pure PBB®<sup>15-17</sup>


• PBB® is an innovative, patented intraocular blue dye component designed to optimise safety profile, with:<sup>15-17</sup>

- ✓ 99% purity<sup>15</sup>
- ✓ Greater ILM tissue selectivity compared to BBG<sup>16,17</sup>
- ✓ Lower retinal distribution versus BBG<sup>16,17</sup>
- ✓ Minimal dye reaching inner retinal layers thus reducing risk of iatrogenic retinal damage<sup>16,17</sup>


*Ex vivo* dye distribution in ILM and retina of porcine eyes<sup>16,17</sup>



BBG: Brilliant Blue G; ERM: Epiretinal membrane; ILM: Internal limiting membrane; PBB®: Pure Benzyl-Brilliant Blue.



# BLutein™ Vitreo Lutein Vitreous Staining Dye



- 2% crystalline Lutein

### Designed to:


- Improve visualisation of the vitreous during surgery

### With the aim of:


- Reducing surgery time
- Minimising risk of complications



DYE300  
1 ml sterile single use glass vials  
(10 units per box)



# BLutein™ Single Lutein Blue ILM Staining Dye



- 1% soluble Lutein
- 0.05% PBB®

### Designed to:

- Provide precise identification of the ILM during surgery
- Improve visualisation of the ILM

### With the aim of:


- Reducing surgery time
- Minimising risk of complications




DYE400  
0.5 ml sterile glass pre-filled syringes  
(5 units per box)



ILM: Internal limiting membrane; PBB®: Pure Benzyl-Brilliant Blue.



# BLutein™ Double Lutein Blue ILM and ERM Staining Dye



- 2% soluble Lutein
- 0.05% PBB®
- 0.15% trypan blue

### Designed to:

- Provide precise identification of the ILM and ERM during surgery
- Improve visualisation during surgery

### With the aim of:


- Reducing surgery time
- Minimising risk of complications



DYE500  
0.5 ml sterile glass pre-filled syringes  
(5 units per box)




ERM: Epiretinal membrane; ILM: Internal limiting membrane; PBB®: Pure Benzyl-Brilliant Blue.





# BLutein™ Phaco Lutein Anterior Capsule Staining Dye




- 1% soluble Lutein
- 0.04% trypan blue

### Designed to:

- Provide optimal anterior capsule staining for cataract surgeries
- Facilitate visualisation of the anterior lens capsule
- Facilitate continuous curvilinear capsulorhexis



DYE200  
0.5 ml sterile glass pre-filled syringes  
(5 units per box)



# BLutein™ Dyes

The first and only complete dyes range for posterior and anterior eye surgeries, with the benefits of lutein, designed to help free you from safety profile concerns

- With patented lutein and PBB® selected for their optimised safety profile and proven efficacy<sup>14-17</sup>



## References

1. Bergamo VC *et al.*, Asia Pac J Ophthalmol (Phila). 2021; 10: 26-38. 2. Rodrigues EB *et al.*, Am J Ophthalmol. 2010; 149: 265-77. 3. Casaroli-Marano RP *et al.*, Curr Eye Res. 2015; 40: 707-718. 4. Morales MC *et al.*, Invest Ophthalmol Vis Sci. 2010; 51: 6018-6029. 5. Buscemi S *et al.*, Nutrients. 2018; 10: 1321. 6. Bernstein PS *et al.*, Prog Retin Eye Res. 2016; 50: 34-66. 7. Lima VC *et al.*, Int J Retin Vit. 2016; 2:19. 8. Beccera MO *et al.*, J Funct Foods. 2020; 66: 103771. 9. Junghans A *et al.*, Arch Biochem Biophys. 2001; 391:160-164. 10. Sasaki M *et al.*, J Nutr Biochem. 2012; 23: 423-429. 11. Bian Q *et al.*, Free Radic Biol Med. 2012; 53: 1298-1307. 12. Sundelin SP *et al.*, Free Radic Biol Med. 2001; 31: 217-225. 13. Kim SR *et al.*, Exp Eye Res. 2006; 82: 828-839. 14. Romano MR *et al.*, Graefes Arch Clin Exp Ophthalmol. 2018; 256: 1573-1580. 15. European Patent EP 3 692 101 B1. 16. Spadaro A *et al.*, Frontiers Pharmacol. 2020; 11: 708. 17. Bucolo C *et al.*, Poster #39-A0113. Presented at ARVO Annual Meeting 2019.

Please read the Instructions for Use (IFU) / Directions for Use (DFU) for important product use and safety information for BLutein™ Dyes.

©2023 Bausch + Lomb Incorporated or its affiliates ©/™ are trademarks of Bausch & Lomb Incorporated or its affiliates. DYE INT BLutein Dyes Brochure\_052023\_01.

