

**BAUSCH + LOMB**

# ANTERIOR

Product Catalog





## TABLE OF CONTENTS

STELLARIS ELITE® 4

PHACO PACKS 8

PHACO NEEDLES AND SLEEVES 16

IRRIGATION/ASPIRATION 20

INSTRUMENTS 30


ACCESSORIES 36



The Stellaris Elite® Anterior Deluxe System reinvents responsiveness and revolutionizes control. This level of efficiency and control provides surgeons with a whole new ophthalmic experience. By automatically reacting to crucial procedural factors as they change, surgeons enjoy exacting control throughout the entire procedure. Adaptive Fluidics precision aspiration control delivers lens removal efficiency by providing exceptional chamber stability by automatically compensating for changes in ocular fluid flow.

Working synergistically with the chamber stability and vacuum efficiency of Adaptive Fluidics, the highly efficient and tightly controlled emulsification of Attune energy management facilitates a responsive surgical environment with real-time control.

LEARN MORE AT  
[bauschsurgical.com](http://bauschsurgical.com)



STELLARIS ELITE®



## Stellaris Elite®

Stellaris Elite® Anterior Deluxe System  
Comes with Stellaris Elite® foot pedal

**Stellaris Elite**  
Vision Enhancement System

Stellaris Elite® Foot Control  
BL2390



---

Stellaris® Attune Ultrasonic Handpiece

BL3170  
with SureLock



Needle and sleeve not included.

---

Stellaris® I/A Handpiece

BL3160  
with SureLock



Needle and sleeve not included.

---

I/A Handpiece


DP9721  
Titanium



---

DigiFlow™ Tubing

D4600A  
Air exchange line



---

Digital Media System (DMS)


BL2300HD  
Video overlay system



Adding convenience and efficiency to procedures, our Phaco packs are configured for your procedural needs. Find the perfect fit for your practice and choose from the selection of components that surgeons rely upon for their quality and dependability.


Easily reordered to keep up to date with your workflow, Phaco packs are a convenient way to add efficiency to your practice. The sterilized, single-use components display the high quality engineering and craftsmanship that surgeons have come to expect from Bausch + Lomb Surgical.

LEARN MORE AT  
[bauschsurgical.com](https://bauschsurgical.com)





# PHACO PACKS



## Basic Vacuum Gravity Packs

### Stellaris® Basic Vacuum Phaco Pack


BL5111

A 1 DP5501 sleeve



### Packs come with the items below.

- 1 Vacuum cassette and tubing
- 2 Test chamber
- 3 Needle wrench




## StableChamber® Vacuum Gravity Packs

### Stellaris® StableChamber® Vacuum Pack

BL5112

**Pack Includes:**

A 2 DP5501 sleeves



### Stellaris® StableChamber® Vacuum Pack with MICST™ Needle

BL5113

**Pack Includes:**

A 1 BL3318S purple straight needle

B 2 BL3118 MICST™ sleeves



### Stellaris® StableChamber® Vacuum Pack with MICST™ Angled Needle

BL5113A

**Pack Includes:**

A 1 BL3318AS purple angled needle

B 2 BL3118 MICST™ sleeves



**Packs come with the items below.**

- 1 Vacuum cassette and StableChamber® tubing
- 2 Screen drape
- 3 Tray cover drape
- 4 Remote control drape
- 5 Test chamber
- 6 Needle wrench
- 7 Plug




## Premium Vacuum Gravity Packs

### Stellaris® Premium Vacuum Phaco Pack

BL5110

**Pack Includes:**

**A** 2 DP5501 sleeves



### Stellaris® Premium Vacuum Phaco Pack with MICSTM Straight Needle

BL5114

**Pack Includes:**

**A** 1 BL3318S purple needle

**B** 2 BL3118 MICSTM sleeves



### Stellaris® Premium Vacuum Phaco Pack with MICSTM Angled Needle

BL5114A

**Pack Includes:**

**A** 1 BL3318AS purple needle

**B** 2 BL3118 MICSTM sleeves



**Packs come with the items below.**

- 1 Vacuum cassette and tubing
- 2 Screen drape
- 3 Tray cover drape
- 4 Remote control drape
- 5 Test chamber
- 6 Needle wrench
- 7 Plug



## Premium Vacuum Gravity Packs

### 2.6 to 2.8mm Vacuum Phaco Pack

BLV2260000

**Pack Includes:**

A 2 DP5511 sleeves



A

### 2.2 to 2.4mm Vacuum Phaco Pack

BLV3220000

**Pack Includes:**

A 2 DP5531 sleeves




A

### 1.8 to 2.0mm Vacuum Pack

BLV4180000

**Pack Includes:**

A 2 BL3118 MICS™ sleeves



A

**Packs come with the items below.**

- 1 Vacuum cassette and tubing
- 2 Screen drape
- 3 Tray cover drape
- 4 Test chamber
- 5 Needle wrench
- 6 Plug





## Basic Adaptive Fluidics™ Packs

### Adaptive Fluidics™ Basic Pack with MICST™ Sleeves

BL5115-4

**Pack Includes:**

- A** 2 BL31118 MICST™ sleeves
- B** Adaptive Fluidics cassette and tubing




### Adaptive Fluidics™ Standard Pack with MICST™ Sleeves and Needle

BL5115-2

**Pack Includes:**

- A** BL3420S MICST™ blue straight needle
- B** 2 BL3124S sleeves
- C** Adaptive Fluidics cassette and tubing with Pulse Reflux



**Above packs also include:**

- 1** Screen drape
- 2** Tray cover drape
- 3** Remote control drape
- 4** Test chamber
- 5** Bottle spike
- 6** Needle wrench
- 7** Plug



## Premium Adaptive Fluidics™ Packs

### Adaptive Fluidics™ Premium Pack

BL5115-1

**Pack Includes:**

- A BL3420S MICS™ blue straight needle
- B 2 BL3124S sleeves



A



B

### Adaptive Fluidics™ Premium MICS™ Pack


BL5115-3

**Pack Includes:**

- A BL3420S MICS™ blue straight needle
- B 2 BL3118 MICS™ sleeves



A



B

**Packs come with the items below.**

- 1 Adaptive Fluidics™ cassette and tubing with StableChamber® and Pulse Reflux
- 2 Screen drape
- 3 Tray cover drape
- 4 Remote control drape
- 5 Test chamber
- 6 Needle wrench
- 7 Bottle spike
- 8 Plug



2

3

4



5



6



7




8

Designed for ultimate precision, manufactured from high grade materials, and available in a wide range of sizes, our selection of Phaco needles and sleeves can suit a variety of needs for many types of practices.

Choose from options including a range of angles, straight needles, MircoFlow®, MicroFlow® Plus, and Cortex with thin tips and various colored threads also available for easy identification in the operating room.

LEARN MORE AT  
[bauschsurgical.com](https://bauschsurgical.com)





## NEEDLES AND SLEEVES



## Phaco Needles

### 30° MICST™ 2.0 Angled Needle

BL3420AS

Single use, sterile, 6/box

BL3420A

Limited use, non-sterile, 6/box

RIS: 2.0mm OD: 0.74-0.95mm ID: 0.5-0.79mm



### 30° MICST™ 2.0 Straight Needle

BL3420S

Single use, sterile, 6/box

BL3420

Limited use, non-sterile, 6/box

RIS: 2.0mm OD: 0.74-0.95mm ID: 0.5-0.79mm



### 30° MICST™ Angled Needle

BL3318AS

Single use, sterile, 6/box

BL3318A

Limited use, non-sterile, 6/box

RIS: 1.8-2.00mm OD: 0.74-0.95mm ID: 0.5-0.79mm



### 30° MICST™ Straight Needle

BL3318S

Single use, sterile, 6/box

BL3318

Limited use, non-sterile, 6/box

RIS: 1.8-2.00mm OD: 0.74-0.95mm ID: 0.5-0.79mm



### 30° Thin Tip 20G Angled Needle

DP8430AS

Single use, sterile, 6/box

DP8430

Limited use, non-sterile, 6/box

RIS: 2.65-2.75mm OD: 0.90mm ID: 0.67mm



### 30° Thin Tip 20G Straight Needle

DP8430S

Single use, sterile, 6/box

DP8430A

Limited use, non-sterile, 6/box

RIS: 2.65-2.75mm OD: 0.90mm ID: 0.67mm



## Phaco Sleeves, Accessories and Kits

Kit with (1) BL3420S MICST™ 2.0mm Blue Straight Needle,  
 (2) DP5531 Green Sleeves, Wrench, Test Chamber and Plug  
 BLP324MX30  
 Sterile, 6 pouches/box



S/U Sleeve and Chamber Kit with 1 Sleeve, 1 Chamber/Pouch  
 DP5501  
 Light blue standard D8109 sleeve and test chamber,  
 Sterile, 12 pouches/box



2.4mm Sleeve  
 BL3124  
 Limited use, non-sterile, 6/box



1.8mm MICST™ Sleeve  
 2 Sleeves/Pouch  
 BL3118  
 Sterile, 12 pouches /box



MicroFlow® 2.2mm Sleeve  
 2 Sleeves/Pouch  
 DP5531  
 Sterile, 12 pouches/box



2.4mm Sleeve  
 2 Sleeves/Pouch  
 BL3124S  
 Sterile, 12 pouches/box



Thin Tip 20G Sleeve  
 2 Sleeves/Pouch  
 DP5511  
 Sterile, 12 pouches/box



Standard Sleeve  
 D8109  
 Non-sterile, 10/box



In order to complement the high levels of efficiency provided by Stellaris Elite<sup>®</sup>, we offer a range of irrigation/aspiration in both single use and reusable options. Our I/A handpieces, including our Capsuleguard and bimanual series, also provide the versatility to work with a variety of cataract platforms in addition to the dual-platform Stellaris Elite<sup>®</sup>.

With a large number of tips, materials and angles to choose from, surgeons from a wide range of practices with varied techniques can appreciate the numerous options available.

LEARN MORE AT  
[bauschsurgical.com](https://www.bauschsurgical.com)



## IRRIGATION/ASPIRATION



## I/A Single Use

### Single-Use I/A MICST™

85795ST

I/A handpiece 45° angled tip, with sleeve for 1.8mm incision, 0.3mm aspiration port, single-use, 12/box



## I/A Single Use CapsuleGuard®

### CapsuleGuard® 45° I/A Handpiece

for all systems

85910S

Silicone tip, single-use, 12/box



### CapsuleGuard® 45° I/A Handpiece for Stellaris®

85910ST

Silicone tip, single-use, 12/box



### CapsuleGuard® Straight I/A Handpiece

for all systems

85912S

Silicone tip, single-use, 12/box



### CapsuleGuard® Straight I/A Handpiece

for Stellaris®

85912ST

Silicone tip, single-use, 12/box



## I/A Single Use CapsuleGuard®

CapsuleGuard® Curved I/A Handpiece  
for all systems

85913S

Silicone tip, single-use, 12/box



CapsuleGuard® MICST™ Straight I/A Handpiece for  
Stellaris®

85914ST

Silicone tip, single-use, 12/box



CapsuleGuard® Curved I/A Handpiece  
for Stellaris®

85913ST

Silicone tip, single-use, 12/box



CapsuleGuard® MICST™ 45° I/A Handpiece  
for all systems

85915S

Silicone tip, single-use, 12/box



CapsuleGuard® MICST™ Straight I/A Handpiece for all  
systems

85914S

Silicone tip, single-use, 12/box



CapsuleGuard® MICST™ 45° I/A Handpiece  
for Stellaris®

85915ST

Silicone tip, single-use, 12/box



## I/A Series Handles

### I/A Handle

BL3160

Stellaris® I/A handpiece with Surelock



Needle and sleeve not included.

### I/A Handpiece

DP9721

I/A handpiece, titanium, luer connectors and knurled grip



## I/A Series Reusable Tips

### Reusable I/A Tip

BL3162

I/A tip 0.4mm straight silicone 20G



### Reusable I/A Tip

BL3165

I/A tip 0.4mm curved metal 17.5G



### Reusable I/A Capsule Polisher Tip

DP9744

I/A capsule polisher 45° Angled, metal



### Reusable I/A Tip

BL3163

I/A tip 0.4mm 45° angled, silicone 20G



### Reusable I/A Tip

BL3166

I/A tip 0.4mm 45° angled, metal 17.5G



### Reusable I/A Capsule Polisher Tip with Silicone Sleeve

DP9745

I/A capsule polisher 45° with silicone sleeve



### Reusable I/A Tip

BL3164

I/A tip 0.4mm straight metal 17.5G





## I/A Series Tips

### Soft Tip I/A Straight

DP9739

Soft I/A tip 6/box Straight, limited reuse



### Soft Tip I/A Curved

DP9740

Soft I/A tip 6/box Curved, limited reuse



### Soft Tip I/A 45° Angled

DP9741

Soft I/A tip 6/box 45°, limited reuse



### I/A Tip Straight

DP9722

Fits handle DP9721, Straight I/A tip, 16G O.D. 0.3mm aspiration port



### I/A Tip Curved

DP9723

Fits handle DP9721, Curved I/A tip, 16G O.D. 0.3mm aspiration port



### I/A Tip 45° Angled

DP9724

Fits handle DP9721, Angled I/A tip, 16G 0.3mm aspiration port



### I/A Tip Modified U Tip

DP9725

Fits handle DP9721, 180° Modified U tip, 16G, 0.3mm open end port



### I/A Tip 90° Angled

DP9730

Fits handle DP9721, 90° Angled aspiration tip, 16G, 0.3mm aspiration port



### I/A Tip 21G Straight

DP9731

Fits handle DP9721, Straight I/A tip, 21G, for use with silicone irrigation sleeve, 0.3mm aspiration port



### I/A Tip 21G 45° Angled

DP9733

Fits handle DP9721, 45° Angled tip, 21G, for use with silicone irrigation sleeves, 0.3mm aspiration port



### I/A Tip 21G Curved

DP9734

Fits handle DP9721, curved I/A tip, 21G, for use with silicone irrigation sleeve, 0.3mm aspiration port



### Adaptor

DP9742

Adaptor for titanium handpiece



### Aspirating Adaptor

DP9726

Male luer aspirating adaptor for use with aspirating cannulas, titanium



### O-Rings

DP9728

Replacement rings for reusable I/A tips



---

### I/A Handpiece Adapter 3" Tubing

MVS1061L

Double lumen tubing, female connector on irrigation line, male on aspiration line



---

### I/A Handpiece 17G Straight

MVS1063

Straight 20mm shaft with tip guard, 0.3mm aspiration port, uses 75mm of double lumen tubing (MVS 1061L)



---

### I/A Handpiece 17G 60°

MVS1063C

Straight 20mm shaft, 2mm angled 60° tip with tip guard, 0.3mm aspiration port, uses 75mm of double lumen tubing (MVS 1061L)



---

### I/A Handpiece 17G Curved

MVS1063S

Curved 20mm shaft, with tip guard, 0.3mm aspiration port, uses 75mm of double lumen tubing (MVS 1061L)



## Bimanual I/A Single Use - 21G Procedure for Stellaris® Systems

### S/U Bimanual 21G Aspiration Handpiece

85780S

Roughened tip, 0.8mm O.D., Curved shaft, 12/box



## Bimanual I/A Single Use - 23G Procedure for Stellaris® Systems

### S/U Bimanual 23G Aspiration Handpiece

85792S

Roughened tip, 0.6mm O.D., Curved shaft, 12/box



### S/U Bimanual 21G Irrigation Handpiece

85787ST

Smooth tip, 0.8mm O.D., Curved shaft, for Stellaris® only, 12/box



### S/U Bimanual 23G Irrigation Handpiece

85807ST

Smooth tip, 0.6mm O.D. Curved shaft, for Stellaris® only, 12/box



### S/U Bimanual 23G Irrigation Handpiece

85807S

Smooth tip, 0.6mm O.D. Curved shaft, 12/box



## Bimanual I/A Single Use - 23G Procedure for other Phaco Systems

---

### S/U Bimanual 23G Aspiration Handpiece

85792S

Roughened tip, 0.6mm O.D, Curved shaft, 12/box



---

### S/U Bimanual 23G Irrigation Handpiece

85807S

SMOOTH tip, 0.6mm O.D, Curved shaft, 12/box



## Reusable Bimanual

### Bimanual 19G Irrigating Handpiece

MVS1082

Designed for small incision phaco, 0.6mm port, Titanium



### Bimanual 23G Aspiration Handpiece

ET4750 AP23

0.4mm roughened tip, Titanium



### Bimanual 18G Aspirating Handpiece

MVS1083

Designed for small incision phaco, .2mm roughened tip, 0.3mm port, Titanium



### Bimanual 21G Irrigating Handpiece

ET4750 I21

Dual 0.4mm irrigating ports, Titanium



### Bimanual 23G Aspiration Handpiece

ET4750 A23

0.4mm smooth tip, Titanium



### Bimanual 23G Aspiration Handpiece

ET4753 A23


0.3mm smooth tip, Titanium



With our large investment in R&D and manufacturing facilities, we not only have precision engineered single-use and reusable instruments for the most performed procedures, but also instruments that surgeons trust for rare and uncommon procedures. That is why so many have come to trust Bausch + Lomb and Storz® day in and day out, no matter what challenges arise.

Our dedication to meeting the evolving needs of the modern surgical practice is evident in our wide range of instruments, designed for reliable performance and trusted efficiency.

LEARN MORE AT  
[bauschsurgical.com](https://bauschsurgical.com)



# INSTRUMENTS



## Single-Use Knives

### Single-Use 15° Stab Guarded Knife

LP7515G

15° sideport guarded, 6/box



### Single-Use 1.5-2.0mm Trapezoid Knife

LP7596

1.5-2.0mm trapezoidal, 6/  
box



### Single-Use 15° Stab Knife

LP7515

15° sideport, 6/box



### Single-Use 2.75mm Angled Slit Knife

LP7559A

2.75mm, angled up, 6/box



### Single-Use 2.4mm Angled Double Bevel Slit Knife

LP7549ADB

2.4mm, angled double  
bevel, 6/box



### Single-Use 1.6-1.8mm Trapezoid Knife

LP7600

1.6-1.8mm trapezoidal,  
6/box



### Single-Use 2.0-2.2mm Trapezoid Knife

LP7601

2.0-2.2mm trapezoidal,  
6/box



### Single-Use 2.85mm Angled Slit Knife

LP7553A

2.85mm, angled up, 6/box



### Single-Use 2.5mm Angled Single Bevel Up Slit Knife

LP7550A

2.5mm, angled single  
bevel up, 6/box



### Single-Use 1.8-2.2-2.8mm Unitrapezoid Knife

LP7602

1.8-2.2-2.8mm uni-trap,  
6/box



### Single-Use 2.65mm Angled Slit Knife

LP7551A

2.65mm, angled up, 6/box



### Single-Use 2.2mm Angled Double Bevel Slit Knife

LP7548ADB

2.2mm, angled double  
bevel, 6/box



### Single-Use 1.1mm Angled Sideport Knife

LP7592

1.1mm, angled, 6/box





## Diamond Knives

Tru-Size™ Angled Trapezoid Diamond Knife 1.8-2.2mm  
E0130 1.8/2.2



Diamond Profile Angled 2.65mm Keratome  
E0130 2.65



Diamond Profile Angled 1.8mm Keratome  
E0130 1.8MM



Diamond Profile Angled 2.7-2.9mm Trapezoid Knife  
E0130 2.75/2.9



Tru-Size™ Angled Trapezoid Diamond Knife 2.0-2.4mm  
E0130 2.0/2.4



Diamond Profile Angled 2.8mm Keratome  
E0130 2.8



Diamond Profile Angled 2.2mm Keratome  
E0130 2.2MM



## Reusable Capsulorhexis Forceps

### Capsulorhexis Forceps Handle

ET8190 H

Capsulorhexis forceps handle titanium (non-rotating)



### Silverstein MICST™ 1.8mm Capsulorhexis Forceps

E2026

Ultra-thin profile, for sub-2.0mm incision, marks at 5.0mm and 6.0mm



### Capsulorhexis Forceps Tip 13mm Shaft

ET1089 T

For use with ET8190 H handle, for 1.4mm incision, curved shaft, 13mm shaft



### Devine Micro Capsulorhexis Forceps

E2020 S

Angled sharp tips, curved shafts, round handle, for 1.8mm incision



### Capsulorhexis Forceps Tip 15mm Shaft

ET1088 TS

For use with ET8190 H handle, for 1.4mm incision, curved shaft, 15mm shaft



### Devgan Micro Capsulorhexis Forceps

ET2012

Ultra fine tip, cross-action forceps, marks at 2.5mm and 5.0mm



### Utrata MICST™ Capsulorhexis Forceps

E2002 M

Angled 13mm shafts, flat handle non-glare finish



---

### Utrata MICST™ Capsulorhexis Forceps

ET1202

13mm shaft, Titanium



---

### Single-Use Capsulorhexis Forceps

#### Single-Use Capsulorhexis Forceps

SUE01-12

Straight jaws, fine sharp tips, 12/box



---

### Wiley Capsulorhexis Forceps

E2002 W

Angled 12mm shafts, iris stop, marks at 2.5mm and 5mm



---

### Single-Use MICST™ Capsulorhexis Forceps

SUE06-12

Fine angled tips for 1.8mm procedures, 12/box



---

### Single-Use Marked Capsulorhexis Forceps

SUE07-12

Straight jaws with marks at 2.75-3.0mm and 5.5-6.0mm, 12/box



---

### Single-Use MICST™ Curved Capsulorhexis Forceps

SUE08-12

MICS™ capsulorhexis forceps curved, 12/box



With the dual-platform Stellaris Elite®, we have whole-eye solutions that are reliable and convenient. Surgeons trust Bausch + Lomb for their entire practice and a variety of procedures.

Our accessories are engineered to give reliable performance and are made with high-quality components. With easy reordering and robust customer service, we ensure that your practice is covered and make it easy for you to grow with Bausch + Lomb Surgical.

LEARN MORE AT  
[bauschsurgical.com](https://bauschsurgical.com)



## ACCESSORIES



## Anterior Vitrectomy Packs - Stellaris®

### Stellaris® Anterior Vitrectomy Pack

BL5612

20G vitrectomy cutter with male adaptor, tubing, infusion sleeve and irrigation cannula

- Need BL3190 Anterior Vitrectomy Adaptor Kit when using this pack with Stellaris®



### 20G Anterior Vitrectomy Cutter

MVS3201-6

20G vitrectomy cutter without irrigation



## Anterior Vitrectomy - Stellaris® PC or Stellaris Elite® Anterior

### Stellaris® Anterior Vitrectomy Pack

BL5612

20G vitrectomy cutter with male adaptor, tubing, infusion sleeve and irrigation cannula 6/box



### 23G Vitrectomy Cutter

BL5623S

23G vitrectomy cutter only with female connector, needs adaptor from BL3190, 6/box



### 20G Vitrectomy Cutter

BL5620

20G vitrectomy cutter only with female connector, needs adaptor from BL3190, 6/box



### 25G Vitrectomy Cutter

BL5625S

25G vitrectomy cutter only with female connector, needs adaptor from BL3190, 6/box



## Anterior Vitrectomy - Stellaris Elite® Posterior

### 20G 7.5K Pneumatic Vitrectomy Cutter

BL5626

20G vitrectomy cutter only with female connector, needs adaptor from BL3190, 6/box



### 25G Stellaris Elite® Cutter

SE5625

25G vitrectomy cutter only with female connector, needs adaptor from BL3190, 6/box



### 23G Stellaris Elite® Cutter

SE5623

23G vitrectomy cutter only with female connector, needs adaptor from BL3190, 6/box



## Anterior Vitrectomy

### Anterior Vitrectomy Adapter Kit

BL3190

21G cannula, drive line, male adaptor, 6/box



### 23G Single Valved Trocar

BL5283

For pars plana anterior vitrectomy



### 25G Single Valved Trocar

BL5285

For pars plana anterior vitrectomy



### S/U Bimanual 22G Irrigation Handpiece

85806ST

Smooth tip, 0.7mm O.D., curved shaft, for Stellaris® only, 12/box



### S/U Bimanual 23G Irrigation Handpiece

85807ST

Smooth tip, 0.6mm O.D., curved shaft, for Stellaris® only, 12/box



### S/U Bimanual 21G Irrigation Handpiece

85787ST

Smooth tip, 0.8mm O.D., curved shaft for Stellaris® only, 12/box






## Maintenance


### Universal Maintenance Kit

CX7120 Kit contains syringe, stopcock and brush  
CX7123 Universal Aspiration Brush  
CX7124 Universal 3-Way Stopcock



### DigiFlow™ Tubing

D4600A  
Air exchange line



### Stellaris® PC Mayo Tray Adaptor

BL4364  
Fits plastic arm for both Stellaris® PC  
and Stellaris ELITE®



### O-Rings

DP9728  
Replacement rings for reusable I/A tips



## Bipolar

### Reusable Bipolar Cord

CX9430  
with lemo connector



### Single-Use Bipolar Pencil 18G Curved

E7918 C  
2-pin connectors, 10/box



### Reusable Bipolar Cord

S2050 B  
2-pin connectors




### Banana Plug Adaptor

CX9404  
reusable



### Single-Use Bipolar Cord

S2070 C  
2-pin connectors, 10/box



### Reusable Bipolar Cord

CX9400



### Single-Use Bipolar Pencil 18G Straight

E7918  
2-pin connectors, 10/box



## Drapes/Covers

---

### Front Panel Drape

DP5300

Clear plastic, 10/box

---

### Mayo Stand Covers

DRA7352

Clear plastic, 24"X40", 50/box

---

### Patient Half Body Incise Drape

DRA1061

Embossed plastic, one fluid pouch, sterile, 20/box

---

### Patient Half-Body Aperture Drape

DRA1059A

Embossed plastic, one fluid pouch, sterile, 20/box

---

### Patient 3/4 Length Aperture Drape

DRA830A

Blue SMS material, two fluid pouches, sterile, 20/box

---

### Patient 3/4 Length Incise Drape

DRA831

Blue SMS material, two fluid pouches, sterile, 20/box



## Standard Back Table Packs

---

### Standard Cataract Procedure Pack

LYN023BSBLB

With full body aperture SMS drape, two fluid pouches and a separate incise film

---

### Standard Vitrectomy Procedure Pack

LYN023VTBLD

With large EENT split drape and separate half body plastic incise drape

---

### Standard ASC Incise Pack

LYN030INSPA

With half-body embossed plastic incise drape, one fluid pouch

---

### Standard ASC Aperture Pack

LYN030APSPA

With half body embossed plastic aperture drape, one fluid pouch



## Trays

### Microsurgical Instrument Tray Small

E7409


Sterilizable plastic tray, capacity 4-7 instruments



### Microsurgical Instrument Tray Double Layer

E7411


Sterilizable plastic tray, capacity 20-30 instruments



### Microsurgical Instrument Tray Large


E7406

Sterilizable plastic tray, capacity 10-20 instruments



### Stellaris® Ultrasonic Handpiece Tray

BL3379



# Index

## Stellaris Elite®

|               |          |
|---------------|----------|
| BL2390.....   | <b>6</b> |
| BL3170.....   | <b>7</b> |
| BL3160.....   | <b>7</b> |
| DP9721.....   | <b>7</b> |
| D4600A.....   | <b>7</b> |
| BL2300HD..... | <b>7</b> |

## Phaco Packs

|                 |           |
|-----------------|-----------|
| BL5111.....     | <b>10</b> |
| BL5112.....     | <b>11</b> |
| BL5113.....     | <b>11</b> |
| BL5113A.....    | <b>11</b> |
| BL5110.....     | <b>12</b> |
| BL5114.....     | <b>12</b> |
| BL5114A.....    | <b>12</b> |
| BLV2260000..... | <b>13</b> |
| BLV3220000..... | <b>13</b> |
| BLV4180000..... | <b>13</b> |
| BL5115-4.....   | <b>14</b> |
| BL5115-2.....   | <b>14</b> |
| BL5115-1.....   | <b>15</b> |
| BL5115-3.....   | <b>15</b> |

## Phaco Needles and Sleeves

|                 |           |
|-----------------|-----------|
| BL3420AS.....   | <b>18</b> |
| BL3420A.....    | <b>18</b> |
| BL3318AS.....   | <b>18</b> |
| BL3318A.....    | <b>18</b> |
| DP8430AS.....   | <b>18</b> |
| DP8430.....     | <b>18</b> |
| BL3420S.....    | <b>18</b> |
| BL3420.....     | <b>18</b> |
| BL3318S.....    | <b>18</b> |
| BL3318.....     | <b>18</b> |
| DP8430S.....    | <b>18</b> |
| DP8430A.....    | <b>18</b> |
| BLP324MX30..... | <b>19</b> |
| DP5501.....     | <b>19</b> |
| BL3124.....     | <b>19</b> |
| BL3124S.....    | <b>19</b> |
| BL3118.....     | <b>19</b> |
| DP5511.....     | <b>19</b> |
| DP5531.....     | <b>19</b> |
| D8109.....      | <b>19</b> |

## Irrigation/Aspiration

|                  |           |
|------------------|-----------|
| 85795ST.....     | <b>22</b> |
| 85910S.....      | <b>22</b> |
| 85910ST.....     | <b>22</b> |
| 85912S.....      | <b>22</b> |
| 85912ST.....     | <b>22</b> |
| 85913S.....      | <b>23</b> |
| 85913ST.....     | <b>23</b> |
| 85914S.....      | <b>23</b> |
| 85914ST.....     | <b>23</b> |
| 85915S.....      | <b>23</b> |
| 85915ST.....     | <b>23</b> |
| BL3160.....      | <b>24</b> |
| DP9721.....      | <b>24</b> |
| BL3162.....      | <b>24</b> |
| BL3163.....      | <b>24</b> |
| BL3164.....      | <b>24</b> |
| BL3165.....      | <b>24</b> |
| BL3166.....      | <b>24</b> |
| DP9744.....      | <b>24</b> |
| DP9745.....      | <b>24</b> |
| DP9739.....      | <b>25</b> |
| DP9740.....      | <b>25</b> |
| DP9741.....      | <b>25</b> |
| DP9722.....      | <b>25</b> |
| DP9723.....      | <b>25</b> |
| DP9724.....      | <b>25</b> |
| DP9725.....      | <b>25</b> |
| DP9730.....      | <b>25</b> |
| DP9731.....      | <b>25</b> |
| DP9733.....      | <b>25</b> |
| DP9734.....      | <b>25</b> |
| DP9742.....      | <b>25</b> |
| DP9726.....      | <b>25</b> |
| DP9728.....      | <b>25</b> |
| MVS 1061L.....   | <b>26</b> |
| MVS 1063.....    | <b>26</b> |
| MVS 1063C.....   | <b>26</b> |
| MVS 1063S.....   | <b>26</b> |
| 85780S.....      | <b>27</b> |
| 85787ST.....     | <b>27</b> |
| 85792S.....      | <b>27</b> |
| 85807ST.....     | <b>27</b> |
| 85807S.....      | <b>27</b> |
| 85792S.....      | <b>28</b> |
| 85807S.....      | <b>28</b> |
| MVS1082.....     | <b>29</b> |
| MVS1083.....     | <b>29</b> |
| ET4750 A23.....  | <b>29</b> |
| ET4750 AP23..... | <b>29</b> |
| ET4750 I21.....  | <b>29</b> |
| ET4753 A23.....  | <b>29</b> |

## Instruments

|                     |           |
|---------------------|-----------|
| LP7515G.....        | <b>32</b> |
| LP7596.....         | <b>32</b> |
| LP7600.....         | <b>32</b> |
| LP7602.....         | <b>32</b> |
| LP7515.....         | <b>32</b> |
| LP7601.....         | <b>32</b> |
| LP7551A.....        | <b>32</b> |
| LP7559A.....        | <b>32</b> |
| LP7553A.....        | <b>32</b> |
| LP7548ADB.....      | <b>32</b> |
| LP7549ADB.....      | <b>32</b> |
| LP7550A.....        | <b>32</b> |
| LP7592.....         | <b>32</b> |
| E0130 1.8/2.2.....  | <b>33</b> |
| E0130 1.8mm.....    | <b>33</b> |
| E0130 2.0/2.4.....  | <b>33</b> |
| E0130 2.2mm.....    | <b>33</b> |
| E0130 2.65.....     | <b>33</b> |
| E0130 2.75/2.9..... | <b>33</b> |
| E0130 2.8.....      | <b>33</b> |
| ET8190 H.....       | <b>34</b> |
| ET1089 T.....       | <b>34</b> |
| ET1088 TS.....      | <b>34</b> |
| E2002 M.....        | <b>34</b> |
| E2026.....          | <b>34</b> |
| E2020 S.....        | <b>34</b> |
| ET2012.....         | <b>34</b> |
| ET1202.....         | <b>35</b> |
| E2002 W.....        | <b>35</b> |
| SUE01-12.....       | <b>35</b> |
| SUE06-12.....       | <b>35</b> |
| SUE07-12.....       | <b>35</b> |
| SUE08-12.....       | <b>35</b> |

## Accessories

|                  |           |
|------------------|-----------|
| BL5612.....      | <b>38</b> |
| MVS3201-6.....   | <b>38</b> |
| BL5612.....      | <b>38</b> |
| BL5620.....      | <b>38</b> |
| BL5623S.....     | <b>38</b> |
| BL5625S.....     | <b>38</b> |
| BL5626.....      | <b>39</b> |
| SE5623.....      | <b>39</b> |
| SE5625.....      | <b>39</b> |
| BL3190.....      | <b>40</b> |
| BL5283.....      | <b>40</b> |
| BL5285.....      | <b>40</b> |
| 85806ST.....     | <b>40</b> |
| 85807ST.....     | <b>40</b> |
| 85787ST.....     | <b>40</b> |
| CX7120.....      | <b>41</b> |
| CX7123.....      | <b>41</b> |
| CX7124.....      | <b>41</b> |
| BL4364.....      | <b>41</b> |
| D4600A.....      | <b>41</b> |
| DP9728.....      | <b>41</b> |
| CX9430.....      | <b>42</b> |
| S2050 B.....     | <b>42</b> |
| S2070 C.....     | <b>42</b> |
| E7918.....       | <b>42</b> |
| E7918 C.....     | <b>42</b> |
| CX9404.....      | <b>42</b> |
| CX9400.....      | <b>42</b> |
| DP5300.....      | <b>43</b> |
| DRA7352.....     | <b>43</b> |
| DRA1061.....     | <b>43</b> |
| DRA1059A.....    | <b>43</b> |
| DRA830A.....     | <b>43</b> |
| DRA831.....      | <b>43</b> |
| LYN023BSBLB..... | <b>44</b> |
| LYN023VTBLD..... | <b>44</b> |
| LYN030INSPA..... | <b>44</b> |
| LYN030APSPA..... | <b>44</b> |
| E7409.....       | <b>45</b> |
| E7406.....       | <b>45</b> |
| E7411.....       | <b>45</b> |
| BL3379.....      | <b>45</b> |

BAUSCH+LOMB SURGICAL SOLUTIONS

# PARTNERING IS OUR PRACTICE.



**eyeTELLIGENCE®**  
Software



**Stellaris Elite®**  
Vision Enhancement System



**Storz®**  
Ophthalmic Instruments



**enVista® Toric**  
Hydrophobic Acrylic  
Toric IOL



**VICTUS®**  
Femtosecond Laser Platform

For more information or to demo these technologies  
contact your Bausch + Lomb representative at 800.338.2020

**BAUSCH + LOMB**

## Important Safety Information for Stellaris Elite® and Accessories

**INDICATIONS:** The Bausch + Lomb Stellaris Elite® vision enhancement system is intended for the emulsification and removal of cataracts, anterior and posterior segment vitrectomy. The system is designed for use in both anterior and posterior segment surgeries. It provides capabilities for phacoemulsification, coaxial and bimanual irrigation/aspiration, bipolar coagulation, vitrectomy, viscous fluid injection/removal and air/fluid exchange operations. The Stellaris Elite® Vision Enhancement System configured with the laser module is additionally intended for retinal photocoagulation and laser trabeculoplasty.

### CONTRAINDICATIONS

#### All Systems:

- Use of accessories not designated by Bausch + Lomb for use with this equipment may result in serious permanent patient injury, adverse surgical outcome, or damage to the equipment.

#### Systems with Laser Module:

- Photocoagulation is not indicated for patients without pigmentation (albino eyes). In addition, Laser Indirect Ophthalmoscope (LIO) is not indicated for cases involving laser photocoagulation within the arcades.

### WARNINGS

#### All Systems:

- Implantable defibrillators present a risk of injury if triggered by a fibrillatory event during intraocular surgery.
- Electromagnetic interaction between the phacoemulsification (phaco) handpiece and an implanted cardiac pacemaker is unlikely but cannot be ruled out.

#### Systems with Laser Module:

- All support personnel who are present during laser treatment must wear appropriate laser protective eyewear.
- DO NOT look directly into the aiming or treatment laser beam.
- Use of unapproved delivery devices may cause inaccurate laser delivery which could result in serious permanent patient injury.

#### When using the VITESSE® handpiece:

- Use only the Entry Site Alignment (ESA) devices provided with the VITESSE® Handpiece Pack (yellow trocar caps). Do not use any ESA with metal components to avoid particulate in the eye.

#### When using the FREEFLOW™ infusion line:

- Do not attempt to administer intraocular gases or viscous fluids using this device.
- The infusion line loop should be created in the horizontal plane.

### GENERAL CAUTIONS FOR SINGLE USE ACCESSORIES:

- Do not re-sterilize or reuse any single use accessories.
- Do not use if package integrity / sterile barrier has been breached or compromised.
- Do not use or attempt to repair damaged single use products.

This is not all you need to know. **Systems with Laser Module: Misuse of the laser system may lead to dangerous situations and severe injuries. All Systems: See the appropriate Operator Manual for detailed directions, proper use, and full risk and safety information. See individual product instructions for use for detailed information on the use of the VITESSE® Handpiece, vitrectomy packs and cutters, and the FREEFLOW™ infusion line.**

**CAUTION:** Federal (U.S.) Law restricts these devices to sale, by or on the order of a physician.

(800) 338-2020

[bauschsurgical.com](http://bauschsurgical.com)

[customerservice@bauschsurgical.com](mailto:customerservice@bauschsurgical.com)



WARNING: Class 4 Laser Product—visible laser radiation.  
Avoid eye or skin exposure to direct or scattered radiation.

**BAUSCH + LOMB**